



Features:

- ASE Anti shock electronic for safety and comfort
- safe and pain free
- clinical tested with 98 % success rate
- no medicaments, no chirurgical interference
- integrate timer
- digital display
- flat and smooth surface allows a straightforward cleaning and sanitize of the device

IO-DRY basic

Healthy re-experience

The non-invasive therapy method against extreme sweating



maximal efficiency

The mains water Iontophoresis (LWI) is one of the effective, safest and cheapest method to treat excessive sweating.



innovative technology

The consequent input of newest technologies offers the development and fabrication of certified medical devices for innovative and natural therapy.



home treatment

Because of the especial easy and safe handling it suits perfectly for therapy at home.

io-dry basic

The MEDlight therapy devices prices through their easy appliance and maximal trustiness by highest efficacy. The classical coflow has the strongest effect through 100 % current supply and therefore it is preferred for foot therapy.

The selection of correct electrode:

Aluminum: If there is a nickel or chrome allergy, we suggest the variant of light metal.

Stainless steel: stainless steel electrodes offers a benefit of easier cleaning and higher cleaning frequency.

Stainless steel is not easily scratched and offers less contact surface for pollution and chalk. Also the danger of discoloration of the electrodes is smaller.

The set includes:

- Iontophoresis device MEDlight io-dry basic
- mains adapter
- suitcase (divisible in two treatment pans)
- connecting cable
- safety-wall wart
- instruction manual



TECHNICAL DATES

control unit

dimensions:	190 x 49 x 137 mm
weight:	0,5 kg
supply voltage:	12 V
max. power input:	500 mA (Sicherung)
power drain:	max. 6 VA
treatment pressure:	4 - 60 V _{DC}
max. power:	35 mA (Sicherung)
treatment power:	0 - 30 mA (5 mA Sicherheitsreserve)
max. output:	225 mW

